



Daily Procedure for Monitoring Temperatures Using Fridge-Tag® 2L



*By performing steps 4 & 5 you would see if the temperature went out of range after you closed the previous day and before midnight.



REFRIGERATOR

IF THE CURRENT OR MINIMUM TEMPERATURE IS COLDER THAN:
2°C (OR 35°F)

IF THE CURRENT OR MAXIMUM TEMPERATURE IS WARMER THAN:
8°C (OR 46°F)
FOR > 60 MINUTES

Immediately

Call the Vaccine Unit at 617-983-6828 under the following conditions:

FREEZER

IF THE CURRENT TEMPERATURE IS WARMER THAN:
-15°C (OR +5°F)
FOR > 60 MINUTES

IF THE MAXIMUM TEMPERATURE IN THE FREEZER IS:
WARMER THAN -10°C (OR 14°F)